

Certificate of Analysis

Product Description:	Human IL-1β Monoclonal Antibody
Catalog Number:	AB1004
Clone Number:	9D6
Lot Number:	20240322
Immunogen:	Recombinant Human Interleukin 1β Protein
Isotype:	IgG2b
Form:	Liquid
Purification Method:	Protein A
Purity:	$\geq 95\%$
Concentration:	2.5mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8$^{\circ}$C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.