

Certificate of Analysis

| | |
|-----------------------------|---|
| Product Description: | Human IFN-α Monoclonal Antibody |
| Catalog Number: | AB1052 |
| Clone Number: | 9G12 |
| Lot Number: | 20231030 |
| Immunogen: | Recombinant Human IFN-α Protein |
| Isotype: | IgG2a |
| Form: | Liquid |
| Purification Method: | Protein G |
| Purity: | $\geq 95\%$ |
| Concentration: | 2.5mg/mL |
| Buffer: | 10mM PBS, pH7.4 |
| Appearance: | Clear or slightly opalescent solution |
| Host: | Mouse |
| Application: | CLIA/EIA(Detection) |
| Storage: | Store at $-20\pm 5^{\circ}\text{C}$ for 2 years. Avoid repeated freeze-thaw cycles. |
| Note: | For diagnostic use only. Not intended for therapeutic use. |