

Certificate of Analysis

| | |
|-----------------------------|---|
| Product Description: | Human IL-4 Monoclonal Antibody |
| Catalog Number: | AB1007 |
| Clone Number: | 8H1 |
| Lot Number: | 20230214 |
| Immunogen: | Recombinant Human Interleukin 4 Protein |
| Isotype: | IgG1 |
| Form: | Liquid |
| Purification Method: | Protein G |
| Purity: | ≥95% (SDS-PAGE) |
| Concentration: | 2.2mg/mL |
| Buffer: | 10mM PBS, pH7.4 |
| Appearance: | Clear or slightly opalescent solution |
| Host: | Mouse |
| Application: | CLIA/EIA(Capture) |
| Storage: | Store at -20±5°C for 2 years. Avoid repeated freeze-thaw cycles. |
| Note: | For diagnostic use only. Not intended for therapeutic use. |