

## Certificate of Analysis

<b>Product Description:</b>	<b>Human IL-4 Monoclonal Antibody</b>
<b>Catalog Number:</b>	<b>AB1008</b>
<b>Clone Number:</b>	<b>18H6</b>
<b>Lot Number:</b>	<b>20230220</b>
<b>Immunogen:</b>	<b>Recombinant Human Interleukin 4 Protein</b>
<b>Isotype:</b>	<b>IgG1</b>
<b>Form:</b>	<b>Liquid</b>
<b>Purification Method:</b>	<b>Protein G</b>
<b>Purity:</b>	<b>≥95% (SDS-PAGE)</b>
<b>Concentration:</b>	<b>2.2mg/mL</b>
<b>Buffer:</b>	<b>10mM PBS</b>
<b>Appearance:</b>	<b>Clear or slightly opalescent solution</b>
<b>Host:</b>	<b>Mouse</b>
<b>Application:</b>	<b>CLIA/EIA(Detection)</b>
<b>Storage:</b>	<b>Store at 2-8°C for 2 years.</b>
<b>Note:</b>	<b>For diagnostic use only. Not intended for therapeutic use.</b>