

Certificate of Analysis

Product Description:	Human IL-8 Monoclonal Antibody
Catalog Number:	AB1013
Clone Number:	6E6
Lot Number:	20240205
Immunogen:	Recombinant Human Interleukin 8 Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95%
Concentration:	1.31mg/mL
Buffer:	10mM PBS, pH 7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Capture)
Storage:	Store at 2 - 8 °C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.