

Certificate of Analysis

Product Description:	Human IL-10 Monoclonal Antibody
Catalog Number:	AB1016
Clone Number:	9C3-1
Lot Number:	20231123
Immunogen:	Recombinant Human Interleukin 10 Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein G
Purity:	≥95%
Concentration:	1.5mg/mL
Buffer:	10mM PBS, pH7.2
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.