

Certificate of Analysis

Product Description:	Human IL-12P70 Monoclonal Antibody
Catalog Number:	AB1017
Clone Number:	5F8
Lot Number:	20240202
Immunogen:	Recombinant Human Interleukin 12P70 Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95%
Concentration:	2.3mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Rat
Application:	CLIA/EIA(Capture)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.