

Certificate of Analysis

Product Description:	Human TNF-α Monoclonal Antibody
Catalog Number:	AB1002
Clone Number:	D08
Lot Number:	20240328
Immunogen:	Recombinant Human TNF-α Protein
Isotype:	IgG2a
Form:	Liquid
Purification Method:	Protein A
Purity:	$\geq 95\%$
Concentration:	2.0mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.