

Certificate of Analysis

Product Description:	Human IFN-α Monoclonal Antibody
Catalog Number:	AB1052
Clone Number:	3G3
Lot Number:	20230427
Immunogen:	Recombinant Human IFN-α Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein G
Purity:	$\geq 95\%$ (SDS-PAGE)
Concentration:	2.1mg/mL
Buffer:	10mM PBS, pH7.2~7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at $-20\pm 5^{\circ}\text{C}$ for 2 years. Avoid repeated freeze-thaw cycles.
Note:	For diagnostic use only. Not intended for therapeutic use.