

Certificate of Analysis

Product Description:	Human IFN-γ Monoclonal Antibody
Catalog Number:	AB1023
Clone Number:	3A9
Lot Number:	20231230
Immunogen:	Recombinant Human IFN-γ Protein
Isotype:	IgG 1
Form:	Liquid
Purification Method:	Protein A
Purity:	$\geq 95\%$
Concentration:	2.4mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Capture)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.