

Certificate of Analysis

Product Description:	Human IFN-γ Monoclonal Antibody
Catalog Number:	AB1024
Clone Number:	1C3
Lot Number:	20231223
Immunogen:	Recombinant Human IFN-γ Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein A
Purity:	$\geq 95\%$
Concentration:	2.3mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA((Detection))
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.