

Certificate of Analysis

Product Description:	Human P-Tau217 Monoclonal Antibody
Catalog Number:	AB1034
Clone Number:	P-Tau217-4
Lot Number:	20240322
Immunogen:	KLH Conjugated Synthetic Peptide
Isotype:	IgG
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95%
Concentration:	2.4mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Rabbit
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.