

## Certificate of Analysis

<b>Product Description:</b>	<b>Human IL-1<math>\beta</math> Recombinant Protein</b>
<b>Catalog Number:</b>	<b>AG1002</b>
<b>Lot Number:</b>	<b>20220919</b>
<b>Form:</b>	<b>Liquid</b>
<b>Purity:</b>	<b><math>\geq 90\%</math> (SDS-PAGE)</b>
<b>Source :</b>	<b>Human</b>
<b>Expression Host:</b>	<b><i>E.coli</i></b>
<b>Concentration:</b>	<b>0.5 mg/ml</b>
<b>Buffer:</b>	<b>10mM PBS, pH 7.4</b>
<b>Application:</b>	<b>According to the needs of the researcher</b>
<b>Storage:</b>	<b>Store at -80°C for long-term. Avoid repeated freeze-thaw.</b>
<b>Note:</b>	<b>For diagnostic use only. Not intended for therapeutic use.</b>