

## Certificate of Analysis

<b>Product Description:</b>	<b>Human IL-2Ra Recombinant Protein</b>
<b>Catalog Number:</b>	<b>AG1003</b>
<b>Lot Number:</b>	<b>20221118</b>
<b>Form:</b>	<b>Liquid</b>
<b>Purity:</b>	<b>≥90% (SDS-PAGE)</b>
<b>Source :</b>	<b>Human</b>
<b>Expression Host:</b>	<b>HEK293 Cells</b>
<b>Concentration:</b>	<b>1.39 mg/mL</b>
<b>Buffer:</b>	<b>10mM PBS, pH 7.4</b>
<b>Application:</b>	<b>According to the needs of the researcher</b>
<b>Storage:</b>	<b>Store at -80°C for long-term. Avoid repeated freeze-thaw.</b>
<b>Note:</b>	<b>For diagnostic use only. Not intended for therapeutic use.</b>