

Certificate of Analysis

Product Description:	Human IFN-α/2b Recombinant Protein
Catalog Number:	AG1011
Lot Number:	20221011
Form:	Liquid
Purity:	$\geq 90\%$ (SDS-PAGE)
Source :	Human
Expression Host:	HEK293 Cells
Concentration:	1mg/ml
Buffer:	10mM PBS, pH 7.4
Application:	According to the needs of the researcher
Storage:	Store at -80°C for long-term. Avoid repeated freeze-thaw.
Note:	For diagnostic use only. Not intended for therapeutic use.