

Certificate of Analysis

Product Description:	Human IL-2 Monoclonal Antibody
Catalog Number:	AB1050
Clone Number:	6D2
Lot Number:	20230824
Immunogen:	Recombinant Human Interleukin 2 Protein
Isotype:	IgG 1
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95%
Concentration:	1.8 mg/mL
Buffer:	10mM PBS, pH7.2
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.