

Certificate of Analysis

Product Description:	Human Aβ42 Monoclonal Antibody
Catalog Number:	AB1041
Clone Number:	B1
Lot Number:	20230330
Immunogen:	Peptide-KLH
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein A
Purity:	$\geq 95\%$ (SDS-PAGE)
Concentration:	2.4mg/mL
Buffer:	10mM PBS
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Capture)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.