

Certificate of Analysis

Product Description:	Human CA15-3 Monoclonal Antibody
Catalog Number:	AB2004
Clone Number:	15B3
Lot Number:	20230713
Immunogen:	Recombinant Human CA15-3 Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95% (SDS-PAGE)
Concentration:	1.0mg/mL
Buffer:	10mM PBS, pH7.0
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at 2-8°C for 2 years.
Note:	For diagnostic use only. Not intended for therapeutic use.