

Certificate of Analysis

Product Description:	Human GPI Monoclonal Antibody
Catalog Number:	AB1025
Clone Number:	14F4
Lot Number:	20240516
Immunogen:	Glucose-6-Phosphateisomerase, GPI
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein G
Purity:	≥95%
Concentration:	1.8mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Capture)
Storage:	Store at -20±5°C for 2 years. Avoid repeated freeze-thaw cycles.
Note:	For diagnostic use only. Not intended for therapeutic use.