

Certificate of Analysis

Product Description:	2019-nCov-N Monoclonal Antibody
Catalog Number:	Covid-19-AB-N019
Clone Number:	N019
Lot Number:	20220303
Immunogen:	2019-nCov-N Protein
Isotype:	IgG2a
Form:	Liquid
Purification Method:	Protein A
Purity:	≥95% (SDS-PAGE)
Concentration:	4.659mg/mL
Buffer:	10mM PBS, pH7.4
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA((Detection)
Storage:	Store at -20±5°C for 2 years. Avoid repeated freeze-thaw cycles.
Note:	For diagnostic use only. Not intended for therapeutic use.