

Certificate of Analysis

Product Description:	2019-nCov-S Monoclonal Antibody
Catalog Number:	AB-S010
Clone Number:	S010
Lot Number:	20210812
Immunogen:	2019-nCov-S Protein
Isotype:	IgG1
Form:	Liquid
Purification Method:	Protein G
Purity:	≥95% (SDS-PAGE)
Concentration:	2mg/mL
Buffer:	10mM PBS, pH7.2
Appearance:	Clear or slightly opalescent solution
Host:	Mouse
Application:	CLIA/EIA(Detection)
Storage:	Store at -20±5°C for 2 years. Avoid repeated freeze-thaw cycles.
Note:	For diagnostic use only. Not intended for therapeutic use.